

Sugden Lab.

### ***Preparation of Electrocompetent cells***

1. Inoculate single colony of DH5 $\alpha$  into 50 ml No Salt LB medium and culture at 37°C O/N.
2. Add 25 ml culture medium into 1 L of pre-incubated No Salt LB medium to 37°C.
3. Grow the cells at 37°C at approximately 200 rpm.
4. Grow cells to an A<sub>600</sub> of 0.6~7.5.

#### **1. Preparation for cell stock**

1. Keep EtOH (>500 ml) at -20°C
2. Chill four 250 ml tubes, microcentrifuge tubes, tips and vortex mixer at 4°C
3. Prepare two ice boxes A(ice/water/NaCl) and B(dry ice/EtOH)
4. Set up the centrifugation machine at 8,000 rpm at 0°C
5. Transfer culture to box A and chill completely for 10 min

#### **2. Cell Stock**

1. Transfer of chilled culture into two 250 ml tubes
2. Centrifuge at 8,000 rpm for 10 min at 0°C and remove the supernatant carefully and vortex pellets mildly.
3. Resuspend all four pellets in a total volume of 200 ml cold 10% glycerol-dH<sub>2</sub>O. Combine all resuspended cells in one 250 ml bottle.
4. Centrifuge at 8,000 rpm for 10 min at 0°C and remove the supernatant carefully and vortex pellets mildly.
5. Resuspend pellets in 150 ml cold 10% glycerol-dH<sub>2</sub>O by vortex mildly.
6. Centrifuge at 8,000 rpm for 10 min at 0°C and remove the supernatant carefully and vortex pellets mildly.
7. Resuspend pellets in 100 ml cold 10% glycerol-dH<sub>2</sub>O by vortex mildly.
8. Centrifuge at 8,000 rpm for 10 min at 0°C and remove the supernatant carefully and vortex pellets mildly.
9. Resuspend pellets in 2 ml cold 10% glycerol-dH<sub>2</sub>O with 1 ml pipette man.
10. Transfer 100  $\mu$ l aliquot into 1.5 ml microcentrifuge tubes
11. Freeze immediately in box B for 5 min
12. Transfer tubes to deep freezer and freeze overnight at -80°C

### ***Electro-Transformation***

(using Bio-Rad Electroporator)

-Note: These are high efficiency cells. 1 pg of pUC19 yielded >1000 colonies when 10% was plated. Be warned and use good judgment when determining how much DNA to  
Last Modified 06/03/05

## Sugden Lab.

transform with.

1. Thaw one tube of cells on ice. Keep cells cold until electroporation.
2. Dispense 50  $\mu$ l homogenous cell prep to a chilled microcentrifuge tube on ice.
3. Add 1  $\mu$ l DNA solution and mix gently.
4. Set Bio-Rad electroporator to 1.8 kV, 200  $\Omega$  resistance, and 25  $\mu$ F capacitance.
5. Transfer the cell and DNA mixture to a chilled 0.2 cm electroporation cuvette sterilized with 70% EtOH..
6. Have 1 ml SOC at r.t. ready in a pipette before electroporation of the sample.
7. Apply one pulse by pressing and holding both buttons simultaneously until the readout flashes 'Ch9' several times and a beep is heard.
8. Add 1 ml SOC immediately and mix by pipetting several times.
9. Transfer the electroporated sample to a microcentrifuge tube.
10. Incubate for 1 h at 37°C with agitation at about 250 rpm. Pre-incubate LB agar plate with appropriate antibiotics at 37°C.
11. Add variable volumes of your transformation mix to LB agar plates and spread the plates.
12. Incubate plates 12~16 h at 37°C