

Sugden Lab.

Transfection by Electroporation

Want the cells to be growing at 5×10^5 .

Want 5×10^6 cells per transfection.

The total volume is 500 ul per transfection

- 1) Wash cuvettes 3 x with 70 % Ethanol and 2 x with PBS.
 - 2) Prepare DNA in eppy tubes.
 - 4) Take the necessary amount of cells and centrifuge at 1000 for 10 min.
 - 5) Resuspend in appropriate amount of media (5×10^6 cells/500ul media).
 - 6) Add 500 uL of cells to each sample of DNA.
 - 7) Then place 500 ul of cells + DNA in the cuvette and electroporate.
 - 8) Placed in a 10 cm dish with 10 mL of media.
- Incubate it for 48hrs.