

Sugden Lab.

GST Protein Prep. Ver.2

- 1) Grow 50ml of culture in LB or TB +antibiotic o/n at 37°C shaker.
- 2) Dilute culture in LB or TB +antibiotic 1:10
- 3) Grow 3hrs at 37°C.
- 4) Induce culture by adding 0.4 mM IPTG final concentration. (For 50 ml final culture add 20 µl of 1M IPTG).
- 5) Grow at 25°C for 1 hr.
- 6) Spin 5 min at 5 K
- 7) Wash pellet with half volume of cold H₂O. (For 50 ml culture use 25 ml H₂O.)
- 8) Wash bacteria again.
- 9) Resuspend in 1 ml of resuspension buffer per 50 ml of culture.

Resuspension Buffer- NETN + protease inhibitors

20 mM Tris pH 8.0
100 mM NaCl
1 mM EDTA
0.5% NP-40 or Triton-X
1 µg/ml Aprotinin
1 µg/ml PMSF
1 µg/ml Benzaminide

Note: Tim has 100x stock of protease inhibitors.

- 10) Sonicate for 2x for 10 seconds in cold room.
- 11) Pellet debris by spinning at 4°C (Easiest way is to put samples in multiple eppendorfs and spin in cold room at max for 2 min.

Sugden Lab.

Binding of Fusion Protein to Beads

- 1) Pipette 400 μ l of GST-Sepharose into Eppendorf
- 2) Spin beads 15 sec @ 8,000 RPM
- 3) Wash beads 2x with 4°C NETN
- 4) Resuspend pellet in 320 μ l of NETN (Final volume will be ~550 μ l) and put into 4 eppendorfs.
- 5) Add 300 μ l of resuspension buffer and 75 μ l of E. Coli. GST Lysate to each tube.
- 6) Incubate for 30 min while rocking in cold room.
- 7) Spin 15 seconds at 8,000 rpm.
- 8) Wash 3x with ~600 μ l of resuspension buffer
- 9) Remove supernatant and resuspend in appropriate volume resuspension buffer. (You can add this directly to SDS load buffer for running on a gel.)