

Sugden Lab.

Immunofluorescent Staining

Note:

- For suspension cells, keep cell concentration less than 1×10^7 /ml.
- Adhesion cells grown on a 22x22mm coverslip can be stained with 150 ul of antibody solution. I usually drop the antibody solution on a piece of parafilm, and put the coverslip onto the drop.
- Be generous with time and buffer spent on blocking and washing.

1. Fix cells with 4% formaldehyde for 15-30 min.
2. Permeabilize cells with 0.05% Triton X-100 at room temperature for at least 10 min.
3. Rinse cells with PBS twice
4. Block cells with 2-5% calf serum in PBS for 30min at room temperature.
5. Stain with primary antibody diluted in 5% CS+PBS
6. Wash with 5% CS+PBS two times
7. Stain with secondary antibody diluted in 5% CS+PBS
8. Wash with 5% CS+PBS three times. (I usually do 5min for each washing)
9. Mount cells on slides with 50% glycerol. For long term storage, you may want to mount cells with solution containing anti-fade reagents (Vector lab sells it).