

Sugden Lab.

Ligation Reactions:

20ul total volume

3:1 to 10:1 insert to vector ratio

Sticky End Ligation: 30 min to 2hrs @ RT

Blunt End Ligation: at least o/n @ RT

Ligation of Sample

1ul T4 DNA Ligase

2ul 10x T4 DNA ligase buffer

(X)ul insert DNA

(Y)ul vector DNA

(17-X-Y)ul dH₂O

Ligation Control

0ul T4 DNA ligase (no ligase control)

2ul 10x T4 DNA ligase buffer

(X)ul insert DNA

(Y)ul vector DNA

(18-X-Y)ul dH₂O

Transformation of Competent E. coli cells:

Thaw aliquot from -80C

Combine 100ul competent E. coli for P3 plasmids and 10ul ligation sample DNA

Incubate on ice 2min

Heat shock @42C for 45 sec to 1 min

Incubate on ice 2min

Add 500ul SOB to sample

Incubate 37C for 40 min to 1 hr w/ shaking

Spin cells @6000rpm, 6min

Pour off supernatant, resuspend cells in residual liquid (~50 to 100 µl)

Plate cells on selective media, incubate 37C for 12-16hrs

Compare colony number on the ligation control plate to the sample plate(s)

Diagnostic of Transformed Cells:

Select several (i.e.8-12) isolatable colonies from the sample plate and inoculate in 5ml LB media with antibiotic selection (per 5ml: 5µl (100mg/ml) ampicillin, 5 µl (50mg/ml) kanamycin)

Incubate 6-8hrs (optimally) with shaking @37C

Conduct Miniprep protocol, conduct diagnostic digest, run on 1% TAE agarose gel to determine if appropriate band sizes are seen

Preparation of a Glycerol Stock:

900ul from a 5ml o/n bacterial culture

100ul 100% glycerol

Store @-80C