## **Picking Clones**

- 1. Transfect 5 ug of the desired plasmid to the desired adherent cell line.
- 2. Transfect in parallel GFP and carrier DNA (total of 5 ug).
- 3. 48hrs after transfection check for transfection efficiency.
- 4. Taking into account the transfection efficiency split cells. 2 plates with  $3x10^4$ ,  $1x10^5$ ,  $3x10^5$  green cells.
- 5. 48hrs after splitting cells, add drug for selection.
- 6. 10 days after adding the selection drug, look for those colonies that have about 100 cells and circle them.
- 7. Let them grow so that there about 1000 cells.
- 8. 2 weeks after adding the selection drug, pick the colonies.
- 9. If you look under the plate, colonies look like gray circles. Under the microscope, they almost fill up the field of view but do not exceed it.

## To pick colonies

- 1. Circle those colonies that will be picked.
- 2. Take off the media and add PBS to wash.
- 3. Take off the PBS.
- 4. Take the tweezers : rinse in 70% ethanol

Shake off the excess

Rinse in PBS

Shake off the excess

- 5. Grab a small piece of whatman paper that is being soaked in trypsin and place it on top of the colony of cells.
- 6. Repeat step 4 and 5 for each colony.
- 7. Look under the microscope. Cells surrounding the whatman paper should look round if they are being ready to be picked.
- 8. Add 2 mL of media to each well in the 6-well plate. Media should contain all the antibiotics.
- 9. Take the tweezer and wash in ethanol and PBS before grabing the whatman paper and dropping it in a well.
- 10. Repeat the wash of the tweezers between colonies of cells that are picked.
- 11. Look at the plate in the air or under the microscope and make sure that most of the cells within the circle have been removed.
- 12. Pipet up and down each well so that you remove most of the cells still on the whatman paper.
- 13. Look at each well to make sure that there are some cells in it.
- 14. Leave the whatman paper in the well.
- 15. Check the wells 4 or 5 days later. Cells should be spreading and growing.
- 16. Different wells will grow at different rates. Expand the ones that grow fast first and let the other ones keep on growing.
- 17. Also keep a back up.